

R-ECO TAPE 1076 / 90% PCR Recycling content / 75 m x 50 mm / Universal packing tape / Dispenser suitable / Solvent free / Made in EU

THE REVOLUTION OF SUSTAINABLE TAPES

- 90% Certified Recycling Content (PCR)
- Excellent adhesion
- Higher yield
- Less waste






Mag. Michael Quester GmbH Anzengruberstr. 3 A-2540 Bad Vöslau www.michaelguester.com office@michaelquester.com



Technical Data sheet

APPLICATION: PACKAGING TAPE - HOT MELT



RECO 1076 is produced with a PCR (post consumer recycled) BOPET film with a recycling content of > 90%

Stocking time

We recommend to use the product within 6 months from the delivery.

Stocking conditions

The product has to be stocked in a dry environment, at a temperature between 15°C and 25°C. The product must be stored protected from dust, humidity and direct sun exposure. It must be stored horizontally and kept in the original packing.

Use

If the storage temperature has been different from what suggested, bring the tape back to that temperature before application.

Before applying the tape, remove any trace of dust or any other dirt on the carton and press the tape down thoroughly for a perfect and time lasting adherence.